

Treatment week	1	2	3	4	5	6	7	8	9	10
Ankle Fracture Weber A										
Ankle Fracture Weber B operativ										
Ankle Fracture Weber C/ Maisonneuve										
Pilon-Tibial										Load increase on the basis of radiological report
Calcaneus operative										Load increase on the basis of radiological report
Calcaneus conservative										Load increase on the basis of radiological report
Metatarsal		VACOpedes 0° sole								
Jones Fraktur										
Hallux		VACOpedes decompression sole								
Arthrodesis Lower Ankle Joint										Load increase on the basis of radiological report
Arthrodesis Upper Ankle Joint										Load increase on the basis of radiological report
Achilles operative/conservative										
Fibular ligaments							VACOankle brace after 3 weeks			

Sole contact  
 Partial load  
 Full load

Storage VACOsplint

Vacoped closed Immobilization with 0°, 15° or 30° Adapter (dorsal and ventral)

Physical therapy

Vacoped closed release movement Upper Ankle Joint with 10°-0°-10° Adapter dorsal