



Stable -  
Comfortable -  
Individual

**VACO<sup>®</sup>pedes**

VACO<sup>®</sup>pedes - the Safe Solution for your feet

**OPED**



**Stable**

- Vacuum cushion and plastic frames effectively stabilize the foot
- Optimal pressure relief in the fore-foot region
- Functional treatment - the talocrural (ankle) joint is freely mobile

**Comfortable**

- Rocker bottom ensures a good walking gait
- Vacuum cushion provides a comfortable foot bed
- Quick and simple to apply thanks to the Velcro closure system
- Terry cloth covers can be washed and replaced



**Individual**

- Vacuum cushion adapts to anatomical situation
- Reversible adjustment, anytime, anywhere, without assistance
- Only two sizes; product can be used on either foot: no left, no right version



**Indications**

- Toe fractures
- Metatarsal fractures
- Hallux valgus/rigidus
- Arthrodeses of the toe joints
- Hammer and claw toes

**Available sizes (Note: VACOpedes can be used on either foot)**

Standard	Shoe sizes	US men 7-13 / women 9+	EUR 41-46
Small	Shoe sizes	US men 5-7 / women 7-9	EUR 36-40

For more information, please visit our Internet site at [www.vacoped.com](http://www.vacoped.com).

For ordering information, please refer to our price list.



**OPED AG** • Sumpfstraße 5 • CH-6312 Steinhausen  
 Fon +41 (0)41/748 53 35 • Fax +41 (0)41/748 53 40  
 E-mail [contact@oped.ch](mailto:contact@oped.ch) • [www.oped.ch](http://www.oped.ch)

**OPED Inc.** • 47 Loring Drive • Framingham, MA 01702 • USA  
 Phone +1/508/879 0762 • Fax +1/508/879 0769 • Toll free +1/866/447 0601  
 Email [contact@opedinc.com](mailto:contact@opedinc.com) • [www.opedinc.com](http://www.opedinc.com)